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# Chemical Breakthrough Times

**MICROFLEX**  
THE MOST TRUSTED NAME IN GLOVES

CHEMICAL	NITRILE								CHLOROPRENE		LATEX	
	Level	Breakthrough Time	Level	Breakthrough Time	Level	Breakthrough Time	Level	Breakthrough Time	Level	Breakthrough Time	Level	Breakthrough Time
<b>Acetaldehyde</b>	0	5 mins			0	0 mins	0	0 mins	0	0 mins		
<b>Acetic Acid (50%)</b>	2	43 mins			3	120 mins	2	35 mins	1	20 mins		
<b>Acetone</b>	0	0 mins	0	0 mins <sup>3</sup>	0	0 mins	0	0 mins	0	3 mins		
<b>Acetonitrile</b>	0	0 mins	0	0 mins <sup>3</sup>	0	0 mins	0	0 mins	0	0 mins		
<b>Acrylamide (40%)</b>	4	120 mins			4	120 mins	4	120 mins	4	120 mins	4	120 mins
<b>Aluminum Nitrate Nonhydrate (10%)</b>	4	120 mins			4	120 mins	4	120 mins	0	0 mins		
<b>Ammonium Hydroxide (29%)</b>	1	16 mins			1	10 mins	0	5 mins	3	107 mins		
<b>Benzene</b>	0	0 mins			0	0 mins	0	0 mins	0	0 mins		
<b>Butyl Acetate</b>	0	0 mins			0	0 mins	0	0 mins	0	0 mins		
<b>Carbon Disulphide</b>			0	0 mins <sup>3</sup>					0	0 mins <sup>3</sup>		
<b>Chloroform (1%)</b>	0	8 mins			0	8 mins	0	8 mins	0	8 mins		
<b>Dichloromethane</b>	0	0 mins	0	0 mins <sup>3</sup>	0	0 mins	0	0 mins	0	0 mins		
<b>Diethylamine</b>			0	0 mins <sup>3</sup>					0	0 mins <sup>3</sup>		
<b>Dimethyl Sulfoxide (DMSO)</b>	0	5 mins			0	1 min	0	0 mins	1	15 mins		
<b>Ethanol (Ethyl Alcohol)</b>	1	10 mins	1	15 mins	1	10 mins	1	10 mins	1	17 mins		
<b>Ethidium Bromide (5%)</b>	4	120 mins	4	120 mins	4	120 mins	4	120 mins	4	120 mins		
<b>Ethyl Acetate</b>	0	0 mins	0	0 mins <sup>3</sup>	0	0 mins	0	0 mins	0	0 mins		
<b>Gluteraldehyde (50%)</b>	4	120 mins			4	120 mins	4	120 mins	4	120 mins		
<b>Hexane</b>	0	3 mins			1	18 mins	0	0 mins	0	0 mins		
<b>Hydrochloric Acid (50%)</b>	4	120 mins			4	120 mins						
<b>Hydrofluoric Acid (40%)</b>	1	26 mins			1	18 mins						
<b>Isopropanol</b>	2	30 mins			2	35 mins	2	36 mins				
<b>Methanol</b>	0	5 mins	0	6 mins <sup>3</sup>	0	5 mins	0	0 mins	1	14 mins <sup>3</sup>		
<b>Methyl Ethyl Ketone (MEK)</b>	0	0 mins			0	0 mins	0	0 mins				
<b>n-Heptane</b>	2	24 mins	0	0 mins <sup>3</sup>	4	120 mins			0	0 mins <sup>3</sup>	0	3 mins
<b>Nitric Acid (50%)</b>	1	16 mins			1	20 mins	1	12 mins				
<b>Periodic Acid (50%)</b>	4	120 mins			4	120 mins	4	120 mins				
<b>Phenol</b>	4	120 mins			4	120 mins	4	120 mins				
<b>Silver Nitrate (0.17N)</b>	4	120 mins			4	120 mins	4	120 mins				
<b>Sodium Hydroxide (50%)</b>	4	120 mins	4	120 mins <sup>1,3</sup>	4	120 mins	4	120 mins	4	120 mins <sup>1,3</sup>	4	120 mins
<b>Sodium Selenate (10%)</b>							4	120 mins				
<b>Sulfuric Acid (50%)</b>							4	120 mins				
<b>Sulfuric Acid (100%)</b>	1	10 mins	4	120 mins <sup>2,3</sup>	1	15 mins	0	7 mins	4	120 mins <sup>2,3</sup>	1	12 mins
<b>Tetrahydrofuran</b>			0	0 mins <sup>3</sup>					0	0 mins <sup>3</sup>		
<b>Toluene</b>	0	1 min	0	0 mins <sup>3</sup>	0	2 mins			0	0 mins <sup>3</sup>	0	1 min

Last updated: 27 May 2009

Third party testing carried out to ASTM F739 methodologies by ARDL (Arkon Rubber Development Lab, USA). The same breakthrough times would result using EN374-3:2003 test methods. This chart is offered as a guide only as the barrier properties of each glove may be affected by differences in material thickness, chemical concentration, temperature, and length of exposure to chemicals. Unless stated, all tests were carried out at a concentration of 99%.

\*The above tests were carried out to a maximum 120 mins. It is unlikely that users will wear a glove continuously for 480 mins (8 hours!) so it was not deemed necessary to test to 480 mins. However, if there is no chemical breakthrough at 120 mins, there is unlikely to be any chemical breakthrough at 480 mins.

## EN374-3:2003 PROTECTION LEVELS

<10 mins Level 0    >10 mins Level 1    >30 mins Level 2    >60 mins Level 3  
>120 mins Level 4    >240 mins Level 5\*    >480 mins Level 6\*

<sup>1</sup> Tested at 40% conc.    <sup>2</sup> Tested at 96% conc.    <sup>3</sup> Tested by TUV Rheinland, Germany to EN374-3 standard

[www.starlabgroup.com](http://www.starlabgroup.com)

## Penetration vs Permeation

**Penetration:** The movement of a chemical and/or micro-organism through the material, pinholes or other imperfections of a glove. **Permeation:** Gloves do not always act as barriers to liquids and can sometimes act as a sponge, soaking up the liquids and holding them against the skin. Therefore, it is necessary to measure the time that the hazardous liquid takes to breakthrough the glove and comes in contact with the skin.

## EN374-2: 2003 Protection Against Liquid Penetration and Micro Organisms

Penetration is tested by an air and/or water leak test and inspected in compliance with the Acceptable Quality Level.

Performance Level	Acceptable Quality Level (AQL)
Level 3	0.65
Level 2	1.5
Level 1	4.0

AQL Level	No. of Gloves Tested	No. of Failures Allowed
0.65 (Level 3)	200	3
1.5 (Level 2)	200	7
4.0 (Level 1)	200	14

For example, on a production run of 100,000 gloves, the table on the left shows how many gloves must be tested, and the maximum number of failures allowed within that sample batch to achieve acceptance.

## EN374-3 Resistance to Permeation by Chemicals

Glove standards can be very confusing. Not least because different manufacturers use different symbols pertaining to different versions of the same standard. Below is a brief outline of the EN374-3 Standard and the symbols used.

EN374-3: 1994	EN374-3: 2003
<p>The 1994 standard did not differentiate between thin disposable gloves designed for incidental exposure to chemical splashes, and thicker gloves intended for immersion. Only four chemicals had to be tested, these could be chosen by the manufacturer, and the outcome of the tests mattered little. Therefore, this pictogram is merely stating that the testing had been done. The <i>i</i> on the pictogram refers users to more detailed test data displayed on box.</p> 	<p>The pictogram with the flask and fumes can now only be shown if a breakthrough time of at least <b>30 mins</b> has been achieved in at least <b>3 of 12 listed chemicals</b> (see table below). The corresponding code letters of those three chemicals should be shown underneath the shield. This test is a significant improvement on the 1994 standard in terms of its value on those seeking protection from chemicals.</p> 
	<p>This pictogram is used if the glove does not pass the 3 chemicals test explained above, but it does fulfil Level 2 of the penetration test EN374-2:2003, ie it has an AQL of at least 1.5</p>  <p>This pictogram highlights the limitations of the chemical barrier properties of thin gauge disposable gloves and emphasises that they're only designed for incidental splashes. The question mark is to remind users that it refers to 'low chemical resistance' or 'water-proof' gloves.</p> <p>However, some manufacturers, like Microflex, undertake testing on a long list of chemicals with excellent breakthrough results, even though the glove does not pass the 3 of 12 for this standard.</p>

### Who's Certifying Your Gloves for PPE Compliance?

**Where is the Notified Body's ID Number on the box for the gloves you're wearing?**

Without a number under the CE mark to identify the Third Party Registrar (eg. 0197 = TUV Rheinland), the glove manufacturer is self-declaring the glove's EN credentials.

Code Letter	Chemical	Class
A	Methanol	Primary alcohol
B	Acetone	Ketone
C	Acetonitrile	Nitrile compound
D	Dichloromethane	Chlorinated paraffin
E	Carbon disulphide	Sulphur containing organic compound
F	Toluene	Aromatic hydrocarbon
G	Diethylamine	Amine
H	Tetrahydrofurane	Heterocyclic and ether compound
I	Ethyl acetate	Ester
J	n-Heptane	Saturated hydrocarbon
K	Sodium Hydroxide (40%)	Inorganic base
L	Sulphuric Acid (96%)	Inorganic mineral acid

### The 12 Chemicals for EN374-3:2003 Compliance

In each chemical class, it appears that the most aggressive chemical has been selected, giving the glove wearer a worse case scenario for chemicals in that particular classification. Therefore, it is quite likely that no thin gauge disposable glove will achieve over 30 mins on three of these chemicals.

**STARLAB GmbH**  
Kurt Fischer Str. 22  
D-22926 Ahrensburg  
Germany  
T: +49 (0)4102 6044 0  
F: +49 (0)4102 6044 20  
info@starlab.de

**STARLAB (UK), Ltd**  
4 Tanners Drive  
Milton Keynes MK14 5NA  
United Kingdom  
T: +44 (0)1908 283800  
F: +44 (0)1908 283802  
info@starlab.co.uk

**STARLAB-FRANCE**  
38 Avenue Henri Barbusse  
92220 Bagneux  
France  
T: +33 (0)1 45 3652 80  
F: +33 (0)1 45 3652 86  
info@starlab-france.com

**STARLAB Srl**  
Via Corelli 36/5  
20134 Milano  
Italy  
T: +39 02 7020 1040  
F: +39 02 7020 1033  
info@starlabitalia.it

**STARLAB India Ltd**  
First Floor, Old No: 322  
New No: 24 Thambu Chetty St  
Chennai 600001, India  
T: +91 44 42130 702 / 703  
F: +91 44 42130 701  
info@starlab.co.in

[www.starlabgroup.com](http://www.starlabgroup.com)